

Introducing the Cochlear™ Wireless Mini Microphone 2+

The World is Every Child's Classroom.

BRING THEM CLOSER TO SOUND WITH COCHLEAR TRUE WIRELESS™ TECHNOLOGY



Hear now. And always



A front row seat to better hearing.


Children with hearing implant solutions may still struggle to hear and understand speech in noisy situations and over distance. Providing them with good, consistent hearing in these challenging circumstances is essential to helping them achieve success in the classroom and beyond. The Cochlear™ Wireless Mini Microphone 2+ can help bring children closer to sound by streaming high-quality audio directly to their sound processor. It is a cost-effective, easy-to-use solution that brings flexibility and freedom to every listening situation.

ULTIMATE PERFORMANCE WITH FULL CONNECTIVITY.

Range up to

25

metres*



Full-day streaming (up to 11 hours)

Directional microphones improve signal-to-noise ratio


Table Mic functionality (omni) when placed horizontally

Full connectivity to audio sources (Line-in), loop systems (built-in Telecoil), and FM systems (Europlug)

Volume lock function, Mute function, Auto-mute function when dropped

SIGNIFICANTLY BETTER SPEECH UNDERSTANDING IN NOISE AND QUIET.


In a classroom setting, the Mini Microphone 2+ provides significantly better speech recognition in both noise and quiet compared to the sound processor alone.¹



* Clear line of sight

TOP CLASS HEARING.

The Mini Microphone 2+ is perfect for a noisy classroom. Thanks to its long range, it can be placed anywhere to help close the gap between the student and speaker. Whether the student is sitting at their desk, or working in groups, the Mini Microphone 2+ automatically adapts to the situation at hand.



When FM is preferred, the Mini Microphone 2+ provides access to FM systems through the use of an FM receiver. The Mini Microphone 2+ may reduce the cost of FM Systems by enabling a single FM receiver to stream bilaterally/bimodally and to multiple users, without limiting bandwidth.

Only one FM receiver required for bilateral/bimodal or multiple users.



FM

Compatible

LEARNING AND LISTENING BEYOND THE CLASSROOM.

As children move from the classroom, to the car, to their best friend's birthday party, the Mini Microphone 2+ will help make sure they don't miss a beat. Perfect for one-to-one and group situations, the small, portable device can easily be used in any environment. Whether it is their parent or teacher's caring voice, or their favourite coach's encouraging words, the child will hear sound clearer and crisper than before.




Hear now. And always

As the global leader in implantable hearing solutions, Cochlear is dedicated to bringing the gift of sound to people with moderate to profound hearing loss. We have helped over 450,000 people of all ages live full and active lives by reconnecting them with family, friends and community.

We aim to give our recipients the best lifelong hearing experience and access to innovative future technologies. For our professional partners, we offer the industry's largest clinical, research and support networks.

That's why more people choose Cochlear than any other hearing implant company.

 **Cochlear Bone Anchored Solutions AB** Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden Tel: +46 31 792 44 00 Fax: +46 31 792 46 95

Regional offices:

Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

Cochlear Americas 13059 E Peakview Avenue, Centennial, CO 80111, USA Tel: +1 303 790 9010 Fax: +1 303 792 9025

Cochlear AG EMEA Headquarters, Peter Merian-Weg 4, 4052 Basel, Switzerland Tel: +41 61 205 8204 Fax: +41 61 205 8205

Cochlear Latinoamérica, S. A. International Business Park Building 3835, Office 403, Panama Pacifico, Panama Tel: +507 830 6220 Fax: +507 830 6218

Local offices:

Cochlear Deutschland GmbH & Co. KG Karl-Wiechert-Allee 76A, 30625 Hannover, Germany Tel: +49 511 542 770 Fax: +49 511 542 7770

Cochlear Europe Ltd 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey, KT15 2HJ, United Kingdom Tel: +44 1932 26 3400 Fax: +44 1932 26 3426

Nihon Cochlear Co Ltd Ochanomizu-Motomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan Tel: +81 3 3817 0241 Fax: +81 3 3817 0245

Cochlear (HK) Limited Unit 1204, 12/F, 297-307 Hennessy Road, Wan Chai, Hong Kong Tel: +852 2530 5773 Fax: +852 2530 5183

Cochlear Medical Device (Beijing) Co Ltd Unit 2208-2212, Germdale Tower B, 91 Jianguo Road, Chaoyang District, Beijing 100022, PR China Tel: +86 10 5909 7800 Fax: +86 10 5909 7900

Cochlear Korea Ltd 1st floor, Cheongwon Building 33, Teheran-ro 8 gil, Gangnam-gu, Seoul, Korea Tel: +82 2 533 4663 Fax: +82 2 533 8408

Cochlear Benelux NV Schalienhoevedreef 20 I, 2800 Mechelen, Belgium Tel: +32 1579 5511 Fax: +32 1579 5500

Cochlear Medical Device Company India PVT Ltd Platina Bldg, Ground Floor, Plot No. C 59, G Block, BKC, Bandra East, Mumbai 400051, India Tel: +91 22 6112 1111 Fax: +91 22 6112 1100

Cochlear Italia SRL Via Larga 33, 40138 Bologna, Italy Tel: +39 051 601 53 11 Fax: +39 051 39 20 62

Cochlear France S.A.S. 135 route de Saint Simon, 31035 Toulouse, France Tel: +33 5 34 63 85 85 (international) Tel: 0805 200 016 (national) Fax: +33 5 34 63 85 80

Cochlear Nordic AB Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden Tel: +46 31 335 14 61 Fax: +46 31 335 14 60

Cochlear Tıbbi Cihazlar ve Sağlık Hizmetleri Ltd. Şti. Cubuklu Mah. Bogazici Cad. Bogazici Plaza No: 6/1 Kavacik, 34805 Beykoz-Istanbul, Turkey Tel: +90 216 538 5900 Fax: +90 216 538 5919

Cochlear Canada Inc 2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada Tel: +1 416 972 5082 Fax: +1 416 972 5083

www.cochlear.com